

For Research Use Only

TNFRSF9/CD137 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50110-1

Capture Antibody Information

Catalog Number:
68754-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2A3B2
Reactivity:
human
GenBank:
BC006196

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 9
Gene ID:
3604

Detection Antibody Information

Catalog Number:
68754-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
1A2A3
Reactivity:
human
GenBank:
BC006196

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 9
Gene ID:
3604

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50110-1 targets TNFRSF9/CD137 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFRSF9 Monoclonal antibody, PBS Only (Capture) 68754-1-PBS (2A3B2). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TNFRSF9 Monoclonal antibody, PBS Only (Detector) 68754-2-PBS (1A2A3). 100 μ g. Concentration 1 mg/mL.

Alternative TNFRSF9/CD137 matched antibody pairs: MP00870-1, MP00870-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50110-1, TNFRSF9/CD137 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68754-1-PBS. Detection antibody: 68754-2-PBS. Standard:Eg0016. Range:1.563-100 ng/mL