

For Research Use Only

EPHA2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50114-1

Capture Antibody Information

Catalog Number:
68759-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2B12B6
Reactivity:
human
GenBank:
BC037166

Conjugate:
Unconjugated
Full name:
EPH receptor A2
Gene ID:
1969

Detection Antibody Information

Catalog Number:
68759-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2C5E5
Reactivity:
human
GenBank:
BC037166

Conjugate:
Unconjugated
Full name:
EPH receptor A2
Gene ID:
1969

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50114-1 targets EPHA2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EPHA2 Monoclonal antibody, PBS Only (Capture) 68759-1-PBS (2B12B6). 100 μ g. Concentration 1 mg/mL.

Detection antibody: EPHA2 Monoclonal antibody, PBS Only (Detector) 68759-2-PBS (2C5E5). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50114-1, EPHA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68759-1-PBS. Detection antibody: 68759-2-PBS. Standard:Eg0146. Range: 0.195-50 ng/mL