

For Research Use Only

NCR3/NKp30 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50119-1

Capture Antibody Information

Catalog Number:
68763-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E7D7

Reactivity:
human

GenBank:
BC052582

Conjugate:
Unconjugated

Full name:
natural cytotoxicity triggering receptor 3

Gene ID:
259197

Detection Antibody Information

Catalog Number:
68763-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F2G2

Reactivity:
human

GenBank:
BC052582

Conjugate:
Unconjugated

Full name:
natural cytotoxicity triggering receptor 3

Gene ID:
259197

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50119-1 targets NCR3/NKp30 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NCR3/NKp30 Monoclonal antibody, PBS Only (Capture) 68763-1-PBS (1E7D7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: NCR3/NKp30 Monoclonal antibody, PBS Only (Detector) 68763-2-PBS (1F2G2). 100 μ g. Concentration 1 mg/mL.

Alternative NCR3/NKp30 matched antibody pairs: MP00833-1, MP00833-2, MP00833-3, MP00833-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50119-1, NCR3/NKp30 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68763-1-PBS. Detection antibody: 68763-2-PBS. Standard: Eg31663. Range: 0.195-50 ng/mL