

For Research Use Only

# DDX43 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50126-1

## Capture Antibody Information

**Catalog Number:**  
68454-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1C3F9  
**Reactivity:**  
Human  
**GenBank:**  
BC066938  
**Immunogen Catalog Number:**  
Ag11484

**Conjugate:**  
Unconjugated  
**Full name:**  
DEAD (Asp-Glu-Ala-Asp) box polypeptide 43  
**Gene ID:**  
55510

## Detection Antibody Information

**Catalog Number:**  
68454-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
2A1H3  
**Reactivity:**  
Human  
**GenBank:**  
BC066938  
**Immunogen Catalog Number:**  
Ag11484

**Conjugate:**  
Unconjugated  
**Full name:**  
DEAD (Asp-Glu-Ala-Asp) box polypeptide 43  
**Gene ID:**  
55510

## Applications

**Tested Applications:**  
Cytometric bead array, Sandwich ELISA

**Range:**  
6.25-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50126-1 targets DDX43 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: DDX43 Monoclonal antibody, PBS Only (Capture) 68454-2-PBS (1C3F9). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDX43 Monoclonal antibody, PBS Only (Detector) 68454-3-PBS (2A1H3). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50126-1, DDX43 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68454-2-PBS. Detection antibody: 68454-3-PBS. Standard: Ag11484. Range: 6.25-100 ng/mL