

For Research Use Only

PRPH2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50146-2

Capture Antibody Information

Catalog Number:
68780-3-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
3D6G10

Reactivity:
Human

GenBank:
BC074720

Immunogen Catalog Number:
Ag12555

Conjugate:
Unconjugated

Full name:
peripherin 2 (retinal degeneration, slow)

Gene ID:
5961

Detection Antibody Information

Catalog Number:
68780-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1F6B3

Reactivity:
Human

GenBank:
BC074720

Immunogen Catalog Number:
Ag12555

Conjugate:
Unconjugated

Full name:
peripherin 2 (retinal degeneration, slow)

Gene ID:
5961

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50146-2 targets PRPH2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: PRPH2 Monoclonal antibody, PBS Only (Capture) 68780-3-PBS (3D6G10). 100 μ g. Concentration 1 mg/ml.

Detection antibody: PRPH2 Monoclonal antibody, PBS Only (Detector) 68780-2-PBS (1F6B3). 100 μ g. Concentration 1 mg/ml.

Alternative PRPH2 matched antibody pairs: MP50146-1

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50146-2, PRPH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68780-3-PBS. Detection antibody: 68780-2-PBS. Standard: Ag12555. Range: 0.781-100 ng/mL.