

For Research Use Only

RNASEH1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50153-1

Capture Antibody Information

Catalog Number:
68787-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3C7F6

Reactivity:
Human

GenBank:
BC002973

Immunogen Catalog Number:
Ag33948

Conjugate:
Unconjugated

Full name:
ribonuclease H1

Gene ID:
246243

Detection Antibody Information

Catalog Number:
68787-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
1F3A1

Reactivity:
Human

GenBank:
BC002973

Immunogen Catalog Number:
Ag7990

Conjugate:
Unconjugated

Full name:
ribonuclease H1

Gene ID:
246243

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50153-1 targets RNASEH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: RNASEH1 Monoclonal antibody, PBS Only (Capture) 68787-1-PBS (3C7F6). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RNASEH1 Monoclonal antibody, PBS Only (Detector) 68787-2-PBS (1F3A1). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50153-1, RNASEH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68787-1-PBS. Detection antibody: 68787-2-PBS. Standard: Ag33948. Range: 0.391-100 ng/mL.