

For Research Use Only

CEP290 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50156-1

Capture Antibody Information

Catalog Number:
67716-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1H6F6

Reactivity:
Human

GenBank:
BC008641

Immunogen Catalog Number:
Ag17457

Conjugate:
Unconjugated

Full name:
centrosomal protein 290kDa

Gene ID:
80184

Detection Antibody Information

Catalog Number:
67716-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2H4G5

Reactivity:
Human

GenBank:
BC008641

Immunogen Catalog Number:
Ag17457

Conjugate:
Unconjugated

Full name:
centrosomal protein 290kDa

Gene ID:
80184

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50156-1 targets CEP290 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: CEP290 Monoclonal antibody, PBS Only (Capture) 67716-2-PBS (1H6F6). 100 μ g. Concentration 1 mg/ml.

Detection antibody: CEP290 Monoclonal antibody, PBS Only (Detector) 67716-3-PBS (2H4G5). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50156-1, CEP290 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67716-2-PBS. Detection antibody: 67716-3-PBS. Standard: Ag17457. Range: 1.563-100 ng/mL.