

For Research Use Only

NME1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50159-1

Capture Antibody Information

Catalog Number:
68792-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1C3E8
Reactivity:
Human
GenBank:
BC018994
Immunogen Catalog Number:
Ag28529

Conjugate:
Unconjugated
Full name:
non-metastatic cells 1, protein (NM23A) expressed in
Gene ID:
4830

Detection Antibody Information

Catalog Number:
68792-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C6C7
Reactivity:
Human
GenBank:
BC018994
Immunogen Catalog Number:
Ag28529

Conjugate:
Unconjugated
Full name:
non-metastatic cells 1, protein (NM23A) expressed in
Gene ID:
4830

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50159-1 targets NME1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: NME1 Monoclonal antibody, PBS Only (Capture) 68792-1-PBS (1C3E8). 100 µg. Concentration 1 mg/ml.

Detection antibody: NME1 Monoclonal antibody, PBS Only (Detector) 68792-2-PBS (2C6C7). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50159-1, NME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68792-1-PBS. Detection antibody: 68792-2-PBS. Standard:Ag28529. Range: 0.391-100 ng/mL.