

CDH11 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50193-2

Capture Antibody Information

Catalog Number:
68818-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3D12F4

Reactivity:
Human

GenBank:
BC013609

Immunogen Catalog Number:
Ag30142

Conjugate:
Unconjugated

Full name:
cadherin 11, type 2, OB-cadherin (osteoblast)

Gene ID:
1009

Detection Antibody Information

Catalog Number:
68818-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E10B5

Reactivity:
Human

GenBank:
BC013609

Immunogen Catalog Number:
Ag30142

Conjugate:
Unconjugated

Full name:
cadherin 11, type 2, OB-cadherin (osteoblast)

Gene ID:
1009

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50193-2 targets CDH11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CDH11 Monoclonal antibody, PBS Only (Capture/Detector) 68818-2-PBS (3D12F4). 100 µg. Concentration 1 mg/ml.

Detection antibody: CDH11 Monoclonal antibody, PBS Only (Detector) 68818-3-PBS (1E10B5). 100 µg. Concentration 1 mg/ml.

Alternative CDH11 matched antibody pairs: MP50193-1

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50193-2, CDH11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68818-2-PBS. Detection antibody: 68818-3-PBS. Standard:Ag30142. Range: 0.195-100 ng/mL.