

For Research Use Only

C16orf63 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50210-1

Capture Antibody Information

Catalog Number:
68832-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F1B7
Reactivity:
Human
GenBank:
BC022321
Immunogen Catalog Number:
Ag24071

Conjugate:
Unconjugated
Full name:
chromosome 16 open reading frame 63
Gene ID:
123811

Detection Antibody Information

Catalog Number:
68832-2-PBS
Host:
Mouse
Isotype:
IgG3
Purification Method:
Protein A Magarose purification

Clone ID:
2C4C6
Reactivity:
Human
GenBank:
BC022321
Immunogen Catalog Number:
Ag24071

Conjugate:
Unconjugated
Full name:
chromosome 16 open reading frame 63
Gene ID:
123811

Applications

Tested Applications:
Sandwich ELISA

Range:
0.049-50 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50210-1 targets C16orf63 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: C16orf63 Monoclonal antibody, PBS Only (Capture) 68832-1-PBS (1F1B7). 100 µg. Concentration 1 mg/ml.

Detection antibody: C16orf63 Monoclonal antibody, PBS Only (Detector) 68832-2-PBS (2C4C6). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptgcn.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50210-1, C16orf63 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68832-1-PBS. Detection antibody: HRP-conjugated 68832-2-PBS. Standard: Ag24071. Range: 0.049-50 ng/mL