

For Research Use Only

BCL6 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50215-1

Capture Antibody Information

Catalog Number:
66340-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
3B1D6

Reactivity:
Human

GenBank:
BC150184

Immunogen Catalog Number:
Ag15519

Conjugate:
Unconjugated

Full name:
B-cell CLL/lymphoma 6

Gene ID:
604

Detection Antibody Information

Catalog Number:
66340-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2H5B4

Reactivity:
Human

GenBank:
BC150184

Immunogen Catalog Number:
Ag15519

Conjugate:
Unconjugated

Full name:
B-cell CLL/lymphoma 6

Gene ID:
604

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50215-1 targets BCL6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCL6 Monoclonal antibody, PBS Only (Capture) 66340-5-PBS (3B1D6). 100 μ g. Concentration 1 mg/ml.

Detection antibody: BCL6 Monoclonal antibody, PBS Only (Capture/Detector) 66340-2-PBS (2H5B4). 100 μ g. Concentration 1 mg/ml.

Alternative BCL6 matched antibody pairs: MP80008-1, MP80008-2, MP80008-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50215-1, BCL6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66340-5-PBS. Detection antibody: 66340-2-PBS. Standard: ag15519. Range: 3.125-100 ng/mL.