

For Research Use Only

# BAP1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50223-1

## Capture Antibody Information

Catalog Number:  
68840-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1G9G1

Reactivity:  
Human

GenBank:  
BC001596

Immunogen Catalog Number:  
Ag27333

Conjugate:  
Unconjugated

Full name:  
BRCA1 associated protein-1 (ubiquitin  
carboxy-terminal hydrolase)

Gene ID:  
8314

## Detection Antibody Information

Catalog Number:  
68840-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2C8E3

Reactivity:  
Human

GenBank:  
BC001596

Immunogen Catalog Number:  
Ag27333

Conjugate:  
Unconjugated

Full name:  
BRCA1 associated protein-1 (ubiquitin  
carboxy-terminal hydrolase)

Gene ID:  
8314

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.391-100 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP50223-1 targets BAP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BAP1 Monoclonal antibody, PBS Only (Capture) 68840-1-PBS (1G9G1). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: BAP1 Monoclonal antibody, PBS Only (Detector) 68840-2-PBS (2C8E3). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative BAP1 matched antibody pairs: MP00357-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50223-1, BAP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68840-1-PBS. Detection antibody: 68840-2-PBS. Standard: Ag27333. Range: 0.391-100 ng/mL.