

TSPAN3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP50243-3

Capture Antibody Information

Catalog Number:
68852-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A10E3

Reactivity:
Human

GenBank:
BC011206

Immunogen Catalog Number:
Ag3689

Conjugate:
Unconjugated

Full name:
tetraspanin 3

Gene ID:
10099

Detection Antibody Information

Catalog Number:
68852-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2A11B3

Reactivity:
Human

GenBank:
BC011206

Immunogen Catalog Number:
Ag3689

Conjugate:
Unconjugated

Full name:
tetraspanin 3

Gene ID:
10099

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-3.125 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50243-3 targets TSPAN3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSPAN3 Monoclonal antibody, PBS Only (Capture) 68852-3-PBS (1A10E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: TSPAN3 Monoclonal antibody, PBS Only (Capture/Detector) 68852-4-PBS (2A11B3). 100 µg. Concentration 1 mg/ml.

Alternative TSPAN3 matched antibody pairs: MP50243-1, MP50243-2, MP50243-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50243-3, TSPAN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68852-3-PBS. Detection antibody: 68852-4-PBS. Standard: Ag3689. Range: 0.098-3.125 ng/mL.