

For Research Use Only

CLCC1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50244-3

Capture Antibody Information

Catalog Number: 68853-3-PBS	Clone ID: 1F9G8	Conjugate: Unconjugated
Host: Mouse	Reactivity: Human	Full name: chloride channel CLIC-like 1
Isotype: IgG1	GenBank: BC002939	Gene ID: 23155
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag24682	

Detection Antibody Information

Catalog Number: 68853-5-PBS	Clone ID: 3D6H8	Conjugate: Unconjugated
Host: Mouse	Reactivity: Human	Full name: chloride channel CLIC-like 1
Isotype: IgG1	GenBank: BC002939	Gene ID: 23155
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag24682	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50244-3 targets CLCC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CLCC1 Monoclonal antibody, PBS Only (Capture) 68853-3-PBS (1F9G8). 100 μ g. Concentration 1 mg/ml.

Detection antibody: CLCC1 Monoclonal antibody, PBS Only (Detector) 68853-5-PBS (3D6H8). 100 μ g. Concentration 1 mg/ml.

Alternative CLCC1 matched antibody pairs: MP50244-1, MP50244-2, MP50244-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50244-3, CLCC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68853-3-PBS. Detection antibody: 68853-5-PBS. Standard:Ag24682. Range: 0.781-100 ng/mL.