

For Research Use Only

# HGS Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50245-3

## Capture Antibody Information

**Catalog Number:**  
67818-4-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
2A9B4  
**Reactivity:**  
Human  
**GenBank:**  
BC003565  
**Immunogen Catalog Number:**  
Ag28610

**Conjugate:**  
Unconjugated  
**Full name:**  
hepatocyte growth factor-regulated tyrosine kinase substrate  
**Gene ID:**  
9146

## Detection Antibody Information

**Catalog Number:**  
67818-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
2C9F6  
**Reactivity:**  
Human  
**GenBank:**  
BC003565  
**Immunogen Catalog Number:**  
Ag28610

**Conjugate:**  
Unconjugated  
**Full name:**  
hepatocyte growth factor-regulated tyrosine kinase substrate  
**Gene ID:**  
9146

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
3.125-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50245-3 targets HGS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HGS Monoclonal antibody, PBS Only (Capture) 67818-4-PBS (2A9B4). 100 µg. Concentration 1 mg/mL.

Detection antibody: HGS Monoclonal antibody, PBS Only (Detector) 67818-3-PBS (2C9F6). 100 µg. Concentration 1 mg/mL.

Alternative HGS matched antibody pairs: MP50245-1, MP50245-2, MP50245-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50245-3, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-4-PBS. Detection antibody: 67818-3-PBS. Standard:Ag28610. Range: 3.125-100 ng/mL.