

For Research Use Only

CHCHD3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50249-1

Capture Antibody Information

Catalog Number:
68259-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2E5B8

Reactivity:
Human

GenBank:
BC014839

Immunogen Catalog Number:
Ag22628

Conjugate:
Unconjugated

Full name:
coiled-coil-helix-coiled-coil-helix domain containing 3

Gene ID:
54927

Detection Antibody Information

Catalog Number:
68259-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2H8H2

Reactivity:
Human

GenBank:
BC014839

Immunogen Catalog Number:
Ag22628

Conjugate:
Unconjugated

Full name:
coiled-coil-helix-coiled-coil-helix domain containing 3

Gene ID:
54927

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50249-1 targets CHCHD3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHCHD3 Monoclonal antibody, PBS Only (Capture) 68259-2-PBS (2E5B8). 100 μ g. Concentration 1 mg/ml.

Detection antibody: CHCHD3 Monoclonal antibody, PBS Only (Detector) 68259-3-PBS (2H8H2). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50249-1, CHCHD3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68259-2-PBS. Detection antibody: 68259-3-PBS. Standard: Ag22628. Range: 1.563-200 ng/mL.