

TIMM8A Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50261-3

Capture Antibody Information

Catalog Number:
68866-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
3B4H2

Reactivity:
Human

GenBank:
BC015093

Immunogen Catalog Number:
Ag33190

Conjugate:
Unconjugated

Full name:
translocase of inner mitochondrial membrane 8 homolog A (yeast)

Gene ID:
1678

Detection Antibody Information

Catalog Number:
68866-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3B8E5

Reactivity:
Human

GenBank:
BC015093

Immunogen Catalog Number:
Ag33190

Conjugate:
Unconjugated

Full name:
translocase of inner mitochondrial membrane 8 homolog A (yeast)

Gene ID:
1678

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50261-3 targets TIMM8A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TIMM8A Monoclonal antibody, PBS Only (Capture) 68866-1-PBS (3B4H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: TIMM8A Monoclonal antibody, PBS Only (Detector) 68866-4-PBS (3B8E5). 100 µg. Concentration 1 mg/ml.

Alternative TIMM8A matched antibody pairs: MP50261-1, MP50261-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50261-3, TIMM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68866-1-PBS. Detection antibody: 68866-4-PBS. Standard: Ag33190. Range: 1.563-100 ng/mL.