

For Research Use Only

# p504S,AMACR Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50263-2

## Capture Antibody Information

<b>Catalog Number:</b> 60240-3-PBS	<b>Clone ID:</b> 5F5D6	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> alpha-methylacyl-CoA racemase
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC009471	<b>Gene ID:</b> 23600
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag8720	

## Detection Antibody Information

<b>Catalog Number:</b> 60240-4-PBS	<b>Clone ID:</b> 1C5B12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> alpha-methylacyl-CoA racemase
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC009471	<b>Gene ID:</b> 23600
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag8720	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 1.563-200 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP50263-2 targets p504S,AMACR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: p504S,AMACR Monoclonal antibody, PBS Only (Capture/Detector) 60240-3-PBS (5F5D6). 100 µg. Concentration 1 mg/ml.

Detection antibody: p504S,AMACR Monoclonal antibody, PBS Only (Detector) 60240-4-PBS (1C5B12). 100 µg. Concentration 1 mg/ml.

Alternative p504S,AMACR matched antibody pairs: MP50263-1, MP50263-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage


**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50263-2, p504S,AMACR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60240-3-PBS. Detection antibody: 60240-4-PBS. Standard:Ag8720. Range: 1.563-200 ng/mL