

STIL Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50266-2

Capture Antibody Information

Catalog Number: 66876-2-PBS	Clone ID: 1G4A3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SCL/TAL1 interrupting locus
Isotype: IgG1	GenBank: BC126223	Gene ID: 6491
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag19124	

Detection Antibody Information

Catalog Number: 66876-4-PBS	Clone ID: 1E4B6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: SCL/TAL1 interrupting locus
Isotype: IgG1	GenBank: BC126223	Gene ID: 6491
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag19124	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50266-2 targets STIL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STIL Monoclonal antibody, PBS Only (Capture) 66876-2-PBS (1G4A3). 100 µg. Concentration 1 mg/ml.

Detection antibody: STIL Monoclonal antibody, PBS Only (Detector) 66876-4-PBS (1E4B6). 100 µg. Concentration 1 mg/ml.

Alternative STIL matched antibody pairs: MP50266-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50266-2, STII Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66876-2-PBS. Detection antibody: 66876-4-PBS. Standard: Ag19124. Range: 1.563-200 ng/mL.