

For Research Use Only

STK32A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50267-3

Capture Antibody Information

Catalog Number: 67528-2-PBS	Clone ID: 1B8C7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: serine/threonine kinase 32A
Isotype: IgG1	GenBank: BC021666	Gene ID: 202374
Purification Method: Protein G purification	Immunogen Catalog Number: Ag10942	

Detection Antibody Information

Catalog Number: 67528-3-PBS	Clone ID: 2B6E5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: serine/threonine kinase 32A
Isotype: IgG2a	GenBank: BC021666	Gene ID: 202374
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag10942	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50267-3 targets STK32A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STK32A Monoclonal antibody, PBS Only (Capture/Detector) 67528-2-PBS (1B8C7). 100 μ g. Concentration 1 mg/ml.

Detection antibody: STK32A Monoclonal antibody, PBS Only (Detector) 67528-3-PBS (2B6E5). 100 μ g. Concentration 1 mg/ml.

Alternative STK32A matched antibody pairs: MP50267-1, MP50267-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50267-3, STK32A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67528-2-PBS. Detection antibody: 67528-3-PBS. Standard: Ag10942. Range: 1.563-200 ng/mL.