

For Research Use Only

MMP3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50277-2

Capture Antibody Information

Catalog Number:
68879-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1H6G6
Reactivity:
human
GenBank:
BC074815

Conjugate:
Unconjugated
Full name:
matrix metallopeptidase 3
(stromelysin 1, progelatinase)
Gene ID:
4314

Detection Antibody Information

Catalog Number:
68879-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1B3G6
Reactivity:
human
GenBank:
BC074815

Conjugate:
Unconjugated
Full name:
matrix metallopeptidase 3
(stromelysin 1, progelatinase)
Gene ID:
4314

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50277-2 targets MMP3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: MMP3 Monoclonal antibody, PBS Only (Capture) 68879-3-PBS (1H6G6). 100 µg. Concentration 1 mg/ml.
Detection antibody: MMP3 Monoclonal antibody, PBS Only (Detector) 68879-4-PBS (1B3G6). 100 µg. Concentration 1 mg/ml.
Alternative MMP3 matched antibody pairs: MP00083-1, MP00083-2, MP00083-3, MP50277-1
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50277-2, MMP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68879-3-PBS. Detection antibody: 68879-4-PBS. Standard:Eg0382. Range: 1.563-100 ng/mL.