

# MSH6 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50291-2

## Capture Antibody Information

<b>Catalog Number:</b> 66172-2-PBS	<b>Clone ID:</b> 5B5A7	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> mutS homolog 6 (E. coli)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC004246	<b>Gene ID:</b> 2956
<b>Purification Method:</b> Protein G purification	<b>Immunogen Catalog Number:</b> Ag12645	

## Detection Antibody Information

<b>Catalog Number:</b> 66172-4-PBS	<b>Clone ID:</b> 2B10B9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Mouse	<b>Reactivity:</b> human	<b>Full name:</b> mutS homolog 6 (E. coli)
<b>Isotype:</b> IgG1	<b>GenBank:</b> BC004246	<b>Gene ID:</b> 2956
<b>Purification Method:</b> Protein G Magarose purification	<b>Immunogen Catalog Number:</b> Ag12645	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 6.25-200 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP50291-2 targets MSH6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MSH6 Monoclonal antibody, PBS Only (Capture) 66172-2-PBS (5B5A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: MSH6 Monoclonal antibody, PBS Only (Detector) 66172-4-PBS (2B10B9). 100 µg. Concentration 1 mg/mL.

Alternative MSH6 matched antibody pairs: MP50291-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50291-2, MSH6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66172-2-PBS. Detection antibody: 66172-4-PBS. Standard: Ag12645. Range: 6.25-200 ng/mL