

Integrin alpha-E Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50294-2

Capture Antibody Information

Catalog Number:
68893-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3H4C10

Reactivity:
human

GenBank:
BC113436

Immunogen Catalog Number:
Ag27410

Conjugate:
Unconjugated

Full name:
integrin, alpha E (antigen CD103,
human mucosal lymphocyte antigen
1; alpha polypeptide)

Gene ID:
3682

Detection Antibody Information

Catalog Number:
68893-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2B10C1

Reactivity:
human

GenBank:
BC113436

Immunogen Catalog Number:
Ag27410

Conjugate:
Unconjugated

Full name:
integrin, alpha E (antigen CD103,
human mucosal lymphocyte antigen
1; alpha polypeptide)

Gene ID:
3682

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP50294-2 targets Integrin alpha-E in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Integrin alpha-E Monoclonal antibody, PBS Only (Capture) 68893-1-PBS (3H4C10). 100 µg. Concentration 1 mg/ml.

Detection antibody: Integrin alpha-E Monoclonal antibody, PBS Only (Detector) 68893-2-PBS (2B10C1). 100 µg. Concentration 1 mg/ml.

Alternative Integrin alpha-E matched antibody pairs: MP50294-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50294-2, Integrin alpha-E Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68893-1-PBS. Detection antibody: 68893-2-PBS. Standard: Ag27410. Range: 0.391-100 ng/mL.