

TRIM11 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50306-1

Capture Antibody Information

Catalog Number: 68904-1-PBS	Clone ID: 2D6D1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tripartite motif-containing 11
Isotype: IgG1	GenBank: BC011629	Gene ID: 81559
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag1308	

Detection Antibody Information

Catalog Number: 68904-2-PBS	Clone ID: 2G5E10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tripartite motif-containing 11
Isotype: IgG1	GenBank: BC011629	Gene ID: 81559
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag1308	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50306-1 targets TRIM11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM11 Monoclonal antibody, PBS Only (Capture) 68904-1-PBS (2D6D1). 100 µg. Concentration 1 mg/ml.

Detection antibody: TRIM11 Monoclonal antibody, PBS Only (Detector) 68904-2-PBS (2G5E10). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50306-1, TRIM11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68904-1-PBS. Detection antibody: 68904-2-PBS. Standard:Ag1308. Range: 0.781-50 ng/mL