

For Research Use Only

ROPN1L Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50320-2

Capture Antibody Information

Catalog Number:
68909-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2H10G2

Reactivity:
human

GenBank:
BC014607

Immunogen Catalog Number:
Ag34015

Conjugate:
Unconjugated

Full name:
ropporin 1-like

Gene ID:
83853

Detection Antibody Information

Catalog Number:
68909-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1F3B9

Reactivity:
human

GenBank:
BC014607

Immunogen Catalog Number:
Ag34015

Conjugate:
Unconjugated

Full name:
ropporin 1-like

Gene ID:
83853

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50320-2 targets ROPN1L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ROPN1L Monoclonal antibody, PBS Only (Capture) 68909-3-PBS (2H10G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: ROPN1L Monoclonal antibody, PBS Only (Detector) 68909-2-PBS (1F3B9). 100 µg. Concentration 1 mg/mL.

Alternative ROPN1L matched antibody pairs: MP50320-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50320-2, ROPN1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68909-3-PBS. Detection antibody: 68909-2-PBS. Standard: Ag34015. Range: 0.781-100 ng/mL.