

For Research Use Only

SEMA7A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50321-2

Capture Antibody Information

Catalog Number:
67397-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H8G1

Reactivity:
human

GenBank:
BC101643

Immunogen Catalog Number:
Ag28804

Conjugate:
Unconjugated

Full name:
semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)

Gene ID:
8482

Detection Antibody Information

Catalog Number:
67397-5-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
3A9D2

Reactivity:
human

GenBank:
BC101643

Immunogen Catalog Number:
Ag28804

Conjugate:
Unconjugated

Full name:
semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)

Gene ID:
8482

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50321-2 targets SEMA7A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SEMA7A Monoclonal antibody, PBS Only (Capture) 67397-4-PBS (1H8G1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: SEMA7A Monoclonal antibody, PBS Only (Detector) 67397-5-PBS (3A9D2). 100 μ g. Concentration 1 mg/ml.

Alternative SEMA7A matched antibody pairs: MP50321-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50321-2, SEMA7A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67397-4-PBS. Detection antibody: 67397-5-PBS. Standard:SEMA7A. Range: 0.391-25 ng/mL