

For Research Use Only

# TSLP Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50325-1

## Capture Antibody Information

Catalog Number:  
68911-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1H1G7

Reactivity:  
human

GenBank:  
BC016720

Immunogen Catalog Number:  
Ag4771

Conjugate:  
Unconjugated

Full name:  
thymic stromal lymphopietin

Gene ID:  
85480

## Detection Antibody Information

Catalog Number:  
68911-2-PBS

Host:  
Mouse

Isotype:  
IgG2a

Purification Method:  
Protein A Magarose purification

Clone ID:  
4G8B11

Reactivity:  
human

GenBank:  
BC016720

Immunogen Catalog Number:  
Ag4771

Conjugate:  
Unconjugated

Full name:  
thymic stromal lymphopietin

Gene ID:  
85480

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50325-1 targets TSLP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSLP Monoclonal antibody, PBS Only (Capture) 68911-1-PBS (1H1G7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: TSLP Monoclonal antibody, PBS Only (Detector) 68911-2-PBS (4G8B11). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative TSLP matched antibody pairs: MP50325-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50325-1, TSLP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68911-1-PBS. Detection antibody: 68911-2-PBS. Standard: Ag4771. Range: 0.098-100 ng/mL.