

For Research Use Only

SLC25A36 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50339-1

Capture Antibody Information

Catalog Number:
67896-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2A4F9

Reactivity:
human

GenBank:
BC014064

Immunogen Catalog Number:
Ag23894

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 36

Gene ID:
55186

Detection Antibody Information

Catalog Number:
67896-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H7D8

Reactivity:
human

GenBank:
BC014064

Immunogen Catalog Number:
Ag23894

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 36

Gene ID:
55186

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50339-1 targets SLC25A36 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A36 Monoclonal antibody, PBS Only (Capture/Detector) 67896-2-PBS (2A4F9). 100 µg. Concentration 1 mg/ml.

Detection antibody: SLC25A36 Monoclonal antibody, PBS Only (Detector) 67896-3-PBS (1H7D8). 100 µg. Concentration 1 mg/ml.

Alternative SLC25A36 matched antibody pairs: MP50339-2, MP50339-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50339-1, SLC25A36 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67896-2-PBS. Detection antibody: 67896-3-PBS. Standard: Ag23894. Range: 0.098-12.5 ng/mL