

For Research Use Only

# SLC25A36 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50339-2

## Capture Antibody Information

Catalog Number:  
67896-4-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1H6C6

Reactivity:  
human

GenBank:  
BC014064

Immunogen Catalog Number:  
Ag23894

Conjugate:  
Unconjugated

Full name:  
solute carrier family 25, member 36

Gene ID:  
55186

## Detection Antibody Information

Catalog Number:  
67896-5-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1C3G1

Reactivity:  
human

GenBank:  
BC014064

Immunogen Catalog Number:  
Ag23894

Conjugate:  
Unconjugated

Full name:  
solute carrier family 25, member 36

Gene ID:  
55186

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50339-2 targets SLC25A36 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A36 Monoclonal antibody, PBS Only (Capture) 67896-4-PBS (1H6C6). 100 µg. Concentration 1 mg/ml.

Detection antibody: SLC25A36 Monoclonal antibody, PBS Only (Detector) 67896-5-PBS (1C3G1). 100 µg. Concentration 1 mg/ml.

Alternative SLC25A36 matched antibody pairs: MP50339-1, MP50339-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50339-2, SLC25A36 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67896-4-PBS. Detection antibody: 67896-5-PBS. Standard: Ag23894. Range: 0.098-50 ng/mL