

PTCH1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50345-2

Capture Antibody Information

Catalog Number: 68918-3-PBS	Clone ID: 4E7B10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: patched homolog 1 (Drosophila)
Isotype: IgG1	GenBank: NM_000264	Gene ID: 5727
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28106	

Detection Antibody Information

Catalog Number: 68918-4-PBS	Clone ID: 3C12D5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: patched homolog 1 (Drosophila)
Isotype: IgG1	GenBank: NM_000264	Gene ID: 5727
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28106	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50345-2 targets PTCH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PTCH1 Monoclonal antibody, PBS Only (Capture) 68918-3-PBS (4E7B10). 100 µg. Concentration 1 mg/ml.

Detection antibody: PTCH1 Monoclonal antibody, PBS Only (Capture/Detector) 68918-4-PBS (3C12D5). 100 µg. Concentration 1 mg/ml.

Alternative PTCH1 matched antibody pairs: MP50345-1, MP50345-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50345-2, PTCH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68918-3-PBS. Detection antibody: 68918-4-PBS. Standard: Ag28106. Range: 0.391-100 ng/mL.