

For Research Use Only

MYC Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50347-2

Capture Antibody Information

Catalog Number:
67447-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G7G1
Reactivity:
human
GenBank:
BC000141
Immunogen Catalog Number:
Ag17519

Conjugate:
Unconjugated
Full name:
v-myc myelocytomatosis viral oncogene homolog (avian)
Gene ID:
4609

Detection Antibody Information

Catalog Number:
67447-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2B5C5
Reactivity:
human
GenBank:
BC000141
Immunogen Catalog Number:
Ag17519

Conjugate:
Unconjugated
Full name:
v-myc myelocytomatosis viral oncogene homolog (avian)
Gene ID:
4609

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50347-2 targets MYC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MYC Monoclonal antibody, PBS Only (Capture) 67447-4-PBS (1G7G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: MYC Monoclonal antibody, PBS Only (Detector) 67447-5-PBS (2B5C5). 100 µg. Concentration 1 mg/mL.

Alternative MYC matched antibody pairs: MP50347-1, MP50347-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50347-2, MYC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67447-4-PBS. Detection antibody: 67447-5-PBS. Standard: Ag17519. Range: 0.781-50 ng/mL