

For Research Use Only

MME,CD10 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50348-2

Capture Antibody Information

Catalog Number:
68921-2-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
3A6B2
Reactivity:
human
GenBank:
BC106070
Immunogen Catalog Number:
Ag12506

Conjugate:
Unconjugated
Full name:
membrane metallo-endopeptidase
Gene ID:
4311

Detection Antibody Information

Catalog Number:
68921-1-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
8D8E6
Reactivity:
human
GenBank:
BC106070
Immunogen Catalog Number:
Ag12506

Conjugate:
Unconjugated
Full name:
membrane metallo-endopeptidase
Gene ID:
4311

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50348-2 targets MME,CD10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MME,CD10 Monoclonal antibody, PBS Only (Capture) 68921-2-PBS (3A6B2). 100 µg. Concentration 1 mg/ml.

Detection antibody: MME,CD10 Monoclonal antibody, PBS Only (Detector) 68921-1-PBS (8D8E6). 100 µg. Concentration 1 mg/ml.

Alternative MME,CD10 matched antibody pairs: MP50348-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50348-2, MME, CD10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68921-2-PBS. Detection antibody: 68921-1-PBS. Standard: Ag12506. Range: 6.25-100 ng/mL