

KLRD1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50356-2

Capture Antibody Information

Catalog Number:
68932-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2B10G12
Reactivity:
human
GenBank:
BC028009
Immunogen Catalog Number:
Ag24477

Conjugate:
Unconjugated
Full name:
killer cell lectin-like receptor subfamily D, member 1
Gene ID:
3824

Detection Antibody Information

Catalog Number:
68932-3-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
3G7D10
Reactivity:
human
GenBank:
BC028009
Immunogen Catalog Number:
Ag24477

Conjugate:
Unconjugated
Full name:
killer cell lectin-like receptor subfamily D, member 1
Gene ID:
3824

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50356-2 targets KLRD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KLRD1 Monoclonal antibody, PBS Only (Capture/Detector) 68932-2-PBS (2B10G12). 100 µg. Concentration 1 mg/ml.

Detection antibody: KLRD1 Monoclonal antibody, PBS Only (Detector) 68932-3-PBS (3G7D10). 100 µg. Concentration 1 mg/ml.

Alternative KLRD1 matched antibody pairs: MP50356-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50356-2, KLRD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68932-2-PBS. Detection antibody: 68932-3-PBS. Standard: Ag24477. Range: 0.781-100 ng/mL.