

For Research Use Only

Synaptophysin Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50359-2

Capture Antibody Information

Catalog Number:
67864-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
4A6B6

Reactivity:
human

GenBank:
BC064550

Immunogen Catalog Number:
Ag11803

Conjugate:
Unconjugated

Full name:
synaptophysin

Gene ID:
6855

Detection Antibody Information

Catalog Number:
67864-5-PBS

Host:
Mouse

Isotype:
IgG3

Purification Method:
Protein A Magarose purification

Clone ID:
2A11B3

Reactivity:
human

GenBank:
BC064550

Immunogen Catalog Number:
Ag11803

Conjugate:
Unconjugated

Full name:
synaptophysin

Gene ID:
6855

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50359-2 targets Synaptophysin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Synaptophysin Monoclonal antibody, PBS Only (Capture) 67864-4-PBS (4A6B6). 100 µg. Concentration 1 mg/ml.

Detection antibody: Synaptophysin Monoclonal antibody, PBS Only (Detector) 67864-5-PBS (2A11B3). 100 µg. Concentration 1 mg/ml.

Alternative Synaptophysin matched antibody pairs: MP00345-1, MP00345-2, MP00345-3, MP00345-4, MP50359-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50359-2, Synaptophysin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67864-4-PBS. Detection antibody: 67864-5-PBS. Standard: Ag11803. Range: 3.125-50 ng/mL.