

For Research Use Only

RIG-1/DDX58 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50370-1

Capture Antibody Information

Catalog Number:
67556-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F5D3
Reactivity:
human
GenBank:
BC132786
Immunogen Catalog Number:
Ag18591

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 58
Gene ID:
23586

Detection Antibody Information

Catalog Number:
67556-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F10G6
Reactivity:
human
GenBank:
BC132786
Immunogen Catalog Number:
Ag18591

Conjugate:
Unconjugated
Full name:
DEAD (Asp-Glu-Ala-Asp) box polypeptide 58
Gene ID:
23586

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50370-1 targets RIG-1/DDX58 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX58 Monoclonal antibody, PBS Only (Capture) 67556-2-PBS (2F5D3). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDX58 Monoclonal antibody, PBS Only (Detector) 67556-3-PBS (2F10G6). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50370-1, RIG-1/DDX58 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67556-2-PBS. Detection antibody: 67556-3-PBS. Standard: Ag18591. Range: 0.391-50 ng/mL