

TPSAB1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50385-1

Capture Antibody Information

Catalog Number: 66174-1-PBS	Clone ID: 6A10G1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tryptase alpha/beta 1
Isotype: IgG1	GenBank: BC028059	Gene ID: 7177
Purification Method: Protein A purification	Immunogen Catalog Number: Ag4167	

Detection Antibody Information

Catalog Number: 66174-2-PBS	Clone ID: 3F11C2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tryptase alpha/beta 1
Isotype: IgG2b	GenBank: BC028059	Gene ID: 7177
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag4167	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50385-1 targets TPSAB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TPSAB1 Monoclonal antibody, PBS Only (Capture) 66174-1-PBS (6A10G1). 100 µg. Concentration 1 mg/ml.

Detection antibody: TPSAB1 Monoclonal antibody, PBS Only (Detector) 66174-2-PBS (3F11C2). 100 µg. Concentration 1 mg/ml.

Alternative TPSAB1 matched antibody pairs: MP00516-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50385-1, TPSAB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66174-1-PBS. Detection antibody: 66174-2-PBS. Standard: Ag4167. Range: 3.125-200 ng/mL.