

For Research Use Only

# ADAM17 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50392-4

## Capture Antibody Information

Catalog Number:  
68725-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1D1G4

Reactivity:  
human

GenBank:  
BC136783

Immunogen Catalog Number:  
Ag34022

Conjugate:  
Unconjugated

Full name:  
ADAM metallopeptidase domain 17

Gene ID:  
6868

## Detection Antibody Information

Catalog Number:  
68725-5-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1E7C6

Reactivity:  
human

GenBank:  
BC136783

Immunogen Catalog Number:  
Ag34022

Conjugate:  
Unconjugated

Full name:  
ADAM metallopeptidase domain 17

Gene ID:  
6868

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50392-4 targets ADAM17 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADAM17 Monoclonal antibody, PBS Only (Capture/Detector) 68725-3-PBS (1D1G4). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: ADAM17 Monoclonal antibody, PBS Only (Detector) 68725-5-PBS (1E7C6). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative ADAM17 matched antibody pairs: MP50392-1, MP50392-2, MP50392-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50392-4, ADAM17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68725-3-PBS. Detection antibody: 68725-5-PBS. Standard:Ag34022. Range: 0.098-100 ng/mL.