

RBMS1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50395-1

Capture Antibody Information

Catalog Number:
68950-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2G4E7

Reactivity:
human

GenBank:
BC018951

Immunogen Catalog Number:
Ag1493

Conjugate:
Unconjugated

Full name:
RNA binding motif, single stranded interacting protein 1

Gene ID:
5937

Detection Antibody Information

Catalog Number:
68950-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2D9G11

Reactivity:
human

GenBank:
BC018951

Immunogen Catalog Number:
Ag1493

Conjugate:
Unconjugated

Full name:
RNA binding motif, single stranded interacting protein 1

Gene ID:
5937

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50395-1 targets RBMS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBMS1 Monoclonal antibody, PBS Only (Capture) 68950-1-PBS (2G4E7). 100 µg. Concentration 1 mg/ml.

Detection antibody: RBMS1 Monoclonal antibody, PBS Only (Detector) 68950-2-PBS (2D9G11). 100 µg. Concentration 1 mg/ml.

Alternative RBMS1 matched antibody pairs: MP00575-1, MP00575-2, MP00575-3, MP50395-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50395-1, RBMS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68950-1-PBS. Detection antibody: 68950-2-PBS. Standard:Ag1493. Range: 1.563-50 ng/mL