

# RBMS1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50395-2

## Capture Antibody Information

Catalog Number:  
68950-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
3E7H1

Reactivity:  
human

GenBank:  
BC018951

Immunogen Catalog Number:  
Ag1493

Conjugate:  
Unconjugated

Full name:  
RNA binding motif, single stranded interacting protein 1

Gene ID:  
5937

## Detection Antibody Information

Catalog Number:  
68950-4-PBS

Host:  
Mouse

Isotype:  
IgG2a

Purification Method:  
Protein A Magarose purification

Clone ID:  
3B6C8

Reactivity:  
human

GenBank:  
BC018951

Immunogen Catalog Number:  
Ag1493

Conjugate:  
Unconjugated

Full name:  
RNA binding motif, single stranded interacting protein 1

Gene ID:  
5937

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50395-2 targets RBMS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBMS1 Monoclonal antibody, PBS Only (Capture) 68950-3-PBS (3E7H1). 100 µg. Concentration 1 mg/ml.

Detection antibody: RBMS1 Monoclonal antibody, PBS Only (Detector) 68950-4-PBS (3B6C8). 100 µg. Concentration 1 mg/ml.

Alternative RBMS1 matched antibody pairs: MP00575-1, MP00575-2, MP00575-3, MP50395-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP50395-2, RBMS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68950-3-PBS. Detection antibody: 68950-4-PBS. Standard: Ag1493. Range: 1.563-100 ng/mL.