

For Research Use Only

CD155/PVR Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50402-1

Capture Antibody Information

Catalog Number: 66913-2-PBS	Clone ID: 1G8B8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: poliovirus receptor
Isotype: IgG2a	GenBank: BC015542	Gene ID: 5817
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag28303	

Detection Antibody Information

Catalog Number: 66913-3-PBS	Clone ID: 1G4A4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: poliovirus receptor
Isotype: IgG1	GenBank: BC015542	Gene ID: 5817
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28303	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50402-1 targets CD155/PVR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD155/PVR Monoclonal antibody, PBS Only (Capture) 66913-2-PBS (1G8B8). 100 µg. Concentration 1 mg/ml.

Detection antibody: CD155/PVR Monoclonal antibody, PBS Only (Detector) 66913-3-PBS (1G4A4). 100 µg. Concentration 1 mg/ml.

Alternative CD155/PVR matched antibody pairs: MP00705-1, MP00705-2, MP00705-3, MP01053-1, MP01053-2, MP01053-3, MP01053-4, MP50402-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50402-1, CD155/PVR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66913-2-PBS. Detection antibody: 66913-3-PBS. Standard: Ag28303. Range: 3.125-100 ng/mL.