

For Research Use Only

TNFRSF14 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50404-2

Capture Antibody Information

Catalog Number:
68953-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D12H9

Reactivity:
human

GenBank:
BC002794

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 14 (herpesvirus entry mediator)

Gene ID:
8764

Detection Antibody Information

Catalog Number:
68953-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A2C11

Reactivity:
human

GenBank:
BC002794

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 14 (herpesvirus entry mediator)

Gene ID:
8764

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50404-2 targets TNFRSF14 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFRSF14 Monoclonal antibody, PBS Only (Capture/Detector) 68953-2-PBS (1D12H9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TNFRSF14 Monoclonal antibody, PBS Only (Detector) 68953-3-PBS (1A2C11). 100 μ g. Concentration 1 mg/mL.

Alternative TNFRSF14 matched antibody pairs: MP00364-1, MP50404-1, MP50404-3, MP50404-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50404-2, TNFRSF 14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68953-2-PBS. Detection antibody: 68953-3-PBS. Standard:Eg0015. Range: 0.195-100 ng/mL