

For Research Use Only

TNFRSF14 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50404-7

Capture Antibody Information

Catalog Number:
68961-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2E11B4

Reactivity:
human

GenBank:
BC002794

Immunogen Catalog Number:
Ag28731

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 14 (herpesvirus entry mediator)

Gene ID:
8764

Detection Antibody Information

Catalog Number:
68961-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2D9B12

Reactivity:
human

GenBank:
BC002794

Immunogen Catalog Number:
Ag28731

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 14 (herpesvirus entry mediator)

Gene ID:
8764

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50404-7 targets TNFRSF14 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFRSF14 Monoclonal antibody, PBS Only (Capture) 68961-3-PBS (2E11B4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: TNFRSF14 Monoclonal antibody, PBS Only (Detector) 68961-4-PBS (2D9B12). 100 μ g. Concentration 1 mg/ml.

Alternative TNFRSF14 matched antibody pairs: MP00364-1, MP50404-1, MP50404-2, MP50404-3, MP50404-4, MP50404-5, MP50404-6

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50404-7, TNFRSF 14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68961-3-PBS. Detection antibody: 68961-4-PBS. Standard: Ag28731. Range: 0.098-12.5 ng/mL