

For Research Use Only

PPIF-Specific Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50406-3

Capture Antibody Information

Catalog Number: 68956-5-PBS	Clone ID: 2B7B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: peptidylprolyl isomerase F
Isotype: IgG2a	GenBank: BC005020	Gene ID: 10105
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag13151	

Detection Antibody Information

Catalog Number: 68956-3-PBS	Clone ID: 3A12B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: peptidylprolyl isomerase F
Isotype: IgG1	GenBank: BC005020	Gene ID: 10105
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag13151	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50406-3 targets PPIF-Specific in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PPIF-Specific Monoclonal antibody, PBS Only (Capture) 68956-5-PBS (2B7B5). 100 μ g. Concentration 1 mg/ml.

Detection antibody: PPIF-Specific Monoclonal antibody, PBS Only (Capture/Detector) 68956-3-PBS (3A12B5). 100 μ g. Concentration 1 mg/ml.

Alternative PPIF-Specific matched antibody pairs: MP50406-1, MP50406-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50406-3, PPIF-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68956-5-PBS. Detection antibody: 68956-3-PBS. Standard: Ag13151. Range: 1.563-100 ng/mL.