

MBTPS2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50410-3

Capture Antibody Information

Catalog Number:
68959-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
3D1A5
Reactivity:
human
GenBank:
BC036465
Immunogen Catalog Number:
Ag4390

Conjugate:
Unconjugated
Full name:
membrane-bound transcription factor
peptidase, site 2
Gene ID:
51360

Detection Antibody Information

Catalog Number:
68959-3-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A Magarose purification

Clone ID:
2B7B4
Reactivity:
human
GenBank:
BC036465
Immunogen Catalog Number:
Ag4390

Conjugate:
Unconjugated
Full name:
membrane-bound transcription factor
peptidase, site 2
Gene ID:
51360

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50410-3 targets MBTPS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MBTPS2 Monoclonal antibody, PBS Only (Capture/Detector) 68959-2-PBS (3D1A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: MBTPS2 Monoclonal antibody, PBS Only (Detector) 68959-3-PBS (2B7B4). 100 µg. Concentration 1 mg/ml.

Alternative MBTPS2 matched antibody pairs: MP50410-1, MP50410-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50410-3, MBTPS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68959-2-PBS. Detection antibody: 68959-3-PBS. Standard:Ag4390. Range: 0.391-100 ng/mL.