

For Research Use Only

SLC25A37 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50411-3

Capture Antibody Information

Catalog Number:
67593-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1C1G1

Reactivity:
human

GenBank:
BC132799

Immunogen Catalog Number:
Ag26216

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 37

Gene ID:
51312

Detection Antibody Information

Catalog Number:
67593-6-PBS

Host:
Mouse

Isotype:
IgG3

Purification Method:
Protein A Magarose purification

Clone ID:
1D11F1

Reactivity:
human

GenBank:
BC132799

Immunogen Catalog Number:
Ag26216

Conjugate:
Unconjugated

Full name:
solute carrier family 25, member 37

Gene ID:
51312

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50411-3 targets SLC25A37 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC25A37 Monoclonal antibody, PBS Only (Capture) 67593-5-PBS (1C1G1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: SLC25A37 Monoclonal antibody, PBS Only (Detector) 67593-6-PBS (1D11F1). 100 μ g. Concentration 1 mg/ml.

Alternative SLC25A37 matched antibody pairs: MP50411-1, MP50411-2, MP50411-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50411-3, SLC25A37 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67593-5-PBS. Detection antibody: 67593-6-PBS. Standard: Ag26216. Range: 6.25-200 ng/mL