

For Research Use Only

# HADHA Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50412-1

## Capture Antibody Information

**Catalog Number:**  
60250-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
3E9B1  
**Reactivity:**  
human, mouse, rat  
**GenBank:**  
BC009235  
**Immunogen Catalog Number:**  
Ag1211

**Conjugate:**  
Unconjugated  
**Full name:**  
hydroxyacyl-Coenzyme A dehydrogenase/3-ketoacyl-Coenzyme A thiolase/enoyl-Coenzyme A hydratase (trifunctional protein), alpha subunit  
**Gene ID:**  
3030

## Detection Antibody Information

**Catalog Number:**  
60250-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
5H5G3  
**Reactivity:**  
human  
**GenBank:**  
BC009235  
**Immunogen Catalog Number:**  
Ag1211

**Conjugate:**  
Unconjugated  
**Full name:**  
hydroxyacyl-Coenzyme A dehydrogenase/3-ketoacyl-Coenzyme A thiolase/enoyl-Coenzyme A hydratase (trifunctional protein), alpha subunit  
**Gene ID:**  
3030

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.563-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50412-1 targets HADHA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: HADHA Monoclonal antibody, PBS Only (Capture) 60250-1-PBS (3E9B1). 100 µg. Concentration 1 mg/mL.  
Detection antibody: HADHA Monoclonal antibody, PBS Only (Detector) 60250-2-PBS (5H5G3). 100 µg. Concentration 1 mg/mL.  
Alternative HADHA matched antibody pairs: MP50412-2  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50412-1, HADHA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60250-1-PBS. Detection antibody: 60250-2-PBS. Standard:Ag1211. Range: 1.563-100 ng/mL.