

DDAH2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50417-2

Capture Antibody Information

Catalog Number: 68154-4-PBS	Clone ID: 1G8G9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: dimethylarginine dimethylaminohydrolase 2
Isotype: IgG1	GenBank: BC001435	Gene ID: 23564
Purification Method: Protein G purification	Immunogen Catalog Number: Ag7127	

Detection Antibody Information

Catalog Number: 68154-1-PBS	Clone ID: 1A10G3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat, rabbit	Full name: dimethylarginine dimethylaminohydrolase 2
Isotype: IgG3	GenBank: BC001435	Gene ID: 23564
Purification Method: Protein A purification	Immunogen Catalog Number: Ag7127	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP50417-2 targets DDAH2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDAH2 Monoclonal antibody, PBS Only (Capture/Detector) 68154-4-PBS (1G8G9). 100 µg. Concentration 1 mg/ml.

Detection antibody: DDAH2 Monoclonal antibody, PBS Only (Detector) 68154-1-PBS (1A10G3). 100 µg. Concentration 1 mg/ml.

Alternative DDAH2 matched antibody pairs: MP50417-1, MP50417-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50417-2, DDAH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68154-4-PBS. Detection antibody: 68154-1-PBS. Standard:Ag7127. Range: 0.391-50 ng/mL