

Integrin alpha-1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50423-2

Capture Antibody Information

Catalog Number:
68211-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2F12E9

Reactivity:
human

GenBank:
BC137121

Immunogen Catalog Number:
Ag30631

Conjugate:
Unconjugated

Full name:
integrin, alpha 1

Gene ID:
3672

Detection Antibody Information

Catalog Number:
68211-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1B8A12

Reactivity:
human

GenBank:
BC137121

Immunogen Catalog Number:
Ag30631

Conjugate:
Unconjugated

Full name:
integrin, alpha 1

Gene ID:
3672

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50423-2 targets Integrin alpha-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Integrin alpha-1 Monoclonal antibody, PBS Only (Capture) 68211-4-PBS (2F12E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: Integrin alpha-1 Monoclonal antibody, PBS Only (Detector) 68211-3-PBS (1B8A12). 100 µg. Concentration 1 mg/ml.

Alternative Integrin alpha-1 matched antibody pairs: MP50423-1, MP50423-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50423-2, Integrin alpha-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68211-4-PBS. Detection antibody: 68211-3-PBS. Standard: Ag30631. Range: 3.125-100 ng/mL.