

For Research Use Only

HMGCR Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50438-3

Capture Antibody Information

Catalog Number:
68975-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2G8G11

Reactivity:
human

GenBank:
BC033692

Immunogen Catalog Number:
Ag4444

Conjugate:
Unconjugated

Full name:
3-hydroxy-3-methylglutaryl-Coenzyme A reductase

Gene ID:
3156

Detection Antibody Information

Catalog Number:
68975-5-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
3B8F2

Reactivity:
human

GenBank:
BC033692

Immunogen Catalog Number:
Ag4444

Conjugate:
Unconjugated

Full name:
3-hydroxy-3-methylglutaryl-Coenzyme A reductase

Gene ID:
3156

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50438-3 targets HMGCR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HMGCR Monoclonal antibody, PBS Only (Capture) 68975-4-PBS (2G8G11). 100 µg. Concentration 1 mg/mL.

Detection antibody: HMGCR Monoclonal antibody, PBS Only (Detector) 68975-5-PBS (3B8F2). 100 µg. Concentration 1 mg/mL.

Alternative HMGCR matched antibody pairs: MP50438-1, MP50438-2

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50438-3, HMGR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68975-4-PBS. Detection antibody: 68975-5-PBS. Standard: Ag4444. Range: 0.781-100 ng/mL.