

MYH1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50449-5

Capture Antibody Information

Catalog Number:
60421-3-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
3E1A3

Reactivity:
human

GenBank:
BC114545

Immunogen Catalog Number:
Ag27992

Conjugate:
Unconjugated

Full name:
myosin, heavy chain 1, skeletal muscle, adult

Gene ID:
4619

Detection Antibody Information

Catalog Number:
60421-4-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
1A4A4

Reactivity:
human

GenBank:
BC114545

Immunogen Catalog Number:
Ag27992

Conjugate:
Unconjugated

Full name:
myosin, heavy chain 1, skeletal muscle, adult

Gene ID:
4619

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50449-5 targets MYH1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MYH1 Monoclonal antibody, PBS Only (Capture/Detector) 60421-3-PBS (3E1A3). 100 μ g. Concentration 1 mg/ml.

Detection antibody: MYH1 Monoclonal antibody, PBS Only (Detector) 60421-4-PBS (1A4A4). 100 μ g. Concentration 1 mg/ml.

Alternative MYH1 matched antibody pairs: MP50449-1, MP50449-2, MP50449-3, MP50449-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50449-5, MYH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60421-3-PBS. Detection antibody: 60421-4-PBS. Standard: Ag27992. Range: 0.098-100 ng/mL.